



Quality Engineering Associates, Inc.

755 Middlesex Turnpike, Unit 3, Billerica, MA 01803 USA

Tel: +1 978-528-2034 • Fax: +1 978-528-2033

Email: info@aea.com • URL: www.aea.com

IMPORTANT BULLETIN FOR PIAS-II USERS!

To ensure smooth operation of older PIAS-IIs running with late-model PCs,

QEA advises these practices:

- Power your device with a self-powered USB 2.0 or 3.0 hub supplying at least 500mA/port.
- Avoid using USB extension cables.