



## Quality Engineering Associates, Inc.

755 Middlesex Turnpike, Unit 3, Billerica, MA 01821 USA

Tel: +1 978-528-2034 Fax: +1 978-528-2033 Email: info@qea.com URL: www.qea.com

### Barcode Reader



The barcode tool decodes the information encoded in barcodes printed on test samples. Key information such as Sample ID, Printer ID, and Test Target ID can be embedded in the barcodes and used to identify not only the sample, printer and target, but also to match the target to a test sequence, or perform other record-keeping functions.

Symbologies supported include:

- 1D Barcodes
  - Codabar
  - Code 128
  - Code 32
  - Code 39
  - Code 93
  - Code 11
  - EAN13
  - EAN8
  - EAN5
  - + 2
  - + 5
  - Interleaved 2 of 5
  - Patch code
  - Postnet
  - UPC-A
  - UPC-B
  - Intelligent Mail
  
- 2D Barcodes
  - DataMatrix
  - PDF417
  - QR Code